



OlyLife Vitality Wand

This device operates using airflow and heating technology. Internal cleanliness is important to maintain performance.

MATERIAL BREAKDOWN



Outer Shell
Plastic



Air Outlet
Heat-resistant plastic



Internal System
Fan-based airflow system

AIR OUTLET CARE

Inspect the air outlet weekly
for dust buildup.

CLEANING METHODS:



- Use a soft dry brush
- Use compressed air with a light spray from a safe distance

NEVER:



- Insert objects into the air outlet
- Spray disinfectant into the vent
- Allow water to enter the air channel



Moisture inside the device may cause internal corrosion.

DAILY CLEANING



- Turn the device off
- Allow the auto-clean cycle to finish
- Unplug the device



CLEANING PROCEDURE:

- Wipe the exterior using a dry microfiber cloth



IF NECESSARY:

- Use a slightly damp cloth with clean water only
- Dry the device immediately



PROFESSIONAL USE GUIDELINES



DO NOT USE THE DEVICE IN:

- Humid bathrooms
- Steam rooms
- Areas near essential oil diffusers



Oil vapor may accumulate on the heating element.



FDA DISCLAIMER: OlyLife products are advanced wellness devices that are safe, non-invasive, and supportive of overall well-being. While PEMF technology has been trusted since the 1970s and certain devices are FDA-cleared for specific wellness uses, OlyLife products are not medical devices and are not intended to diagnose, treat, cure, or prevent any disease.

