


OlyLife SK-M3 Shaken Massager

This device combines fabric components with RF technology, heat, and compression systems.



MATERIAL BREAKDOWN



Main Unit

Plastic housing with internal RF components



Fabric Cover

Removable and washable



Extension Straps

Velcro textile straps



FABRIC COVER CLEANING

RECOMMENDED METHOD:



Hand wash only



Use mild detergent



Water temperature below 30°C



Air dry flat

DO NOT:



Machine dry



Use bleach



Use fabric softener



Iron the fabric



Excessive heat may warp the compression elasticity.



DEVICE BODY CLEANING



AFTER EACH USE:

- Turn off the device
- Wipe the surface using a dry microfiber cloth



IF SWEAT IS PRESENT:

- Use a slightly damp cloth
- Dry immediately



NEVER:

- Spray cleaner directly onto the device
- Submerge the device in water
- Use alcohol on the heating area



STRAP CLEANING

- Wipe straps with a neutral cleaner
- Do not soak the straps
- Do not machine wash



Machine washing can degrade the Velcro quickly.



DEVICE LONGEVITY GUIDELINES

AFTER EACH SESSION:



Allow the device to rest for 20 minutes



Do not recharge while the device is extremely hot



Fully recharge the battery if storing for long periods



DISCLAIMER: OlyLife products are advanced wellness devices that are safe, non-invasive, and supportive of overall well-being. While PEMF technology has been trusted since the 1970s and certain devices are FDA-cleared for specific wellness uses, OlyLife products are not medical devices and are not intended to diagnose, treat, cure, or prevent any disease.

